

# BVI

## Corneal Protection

### A Sterile Single Dose Viscous Transparent Solution



## OP'COVER

[bvimedical.com](http://bvimedical.com)

## OP'COVER Introduction

OP'COVER is intended to temporarily cover the patient's cornea with a protective film during an eye examination with a lens or during eye surgery. OP'COVER avoids dehydration or abrasion of the cornea while maintaining its optical clarity.

### DESCRIPTION:

OP'COVER is a sterile, single-use medical device. OP'COVER comes in a single 2 ml dose of sterile, isotonic, viscous, transparent solution. One box contains 10 single doses. The individual peelable pouch guarantees the external sterility of the single dose.

2.2% m/m of hypromellose in phosphate buffer.  
Does not contain natural latex or preservatives.

## OP'COVER for Multiple Purposes

Corneal protection film for:

- Eye surgeries including laser vision correction
- Examinations with diagnostic lenses

## OP'COVER's Advantages



### FOR THE SURGEON

- **SAVE YOUR TIME**  
one single drop lasts for 20 minutes<sup>1</sup>
- **STAY FOCUSED**  
on your surgical gesture
- **VISUALIZE BETTER**  
thanks to its «magnifying effect»
- **PATIENT SATISFACTION**



### FOR THE PATIENT

- **HYDRATION**  
a drop laid onto the cornea protects the eye from dehydration for 20 minutes<sup>1</sup>
- **REDUCING DISCOMFORT DUE TO ABRASION**  
when surgeon uses lenses
- **IMPROVE POST OPERATIVE COMFORT**  
no itching or irritation reported since 1<sup>st</sup> placement on the market<sup>2</sup>




### FOR CLINIC STAFF

- **EASING THE MANAGEMENT**  
one sterile single-dose for one patient; packaging may stay in the OR
- **TIME SAVING**  
for the sterile nurse
- **TIME GAINING**  
for the entire OR team

## Efficiency

A single drop of OP'COVER versus 10 applications of Balanced Saline Solution.<sup>1</sup>



1. Chen, Y. A., Hirschall, N., & Findl, O. (2011). Comparison of corneal wetting properties of viscous eye lubricant and balanced salt solution to maintain optical clarity during cataract surgery. *Journal of Cataract & Refractive Surgery*, 37(10), 1806-1808.

2. ARCAD\_PMS-PMCF\_11 Version 3 of the 30/09/2022.

## OP'COVER Packaging Features




### HYGROMETRY PRESERVATION:

Resealable bag made of a multi-layer material to preserve single-dose from dehydration

### ASEPSIS PROTECTION:

- Individual pouch guarantees the external sterility of each single-use tube
- Rounded container body to avoid contact of the tip with sterile field



### EASY SQUEEZE:

2 ml plastic tube needing minimal pressure

Qty: 10/Bx

Description	SKU #
OP' Cover	401

## OP'COVER at a Glance



## Safety First for Our Product Line



### ARCIOLANE

Fractionated and Purified  
Silicone Oil




### ARCALINE - ARCOTANE

Heavy Liquids for  
Ophthalmic Surgery



### MONOBLUE

Ultra-purified Staining  
Solutions



### ARCEOLE\*

Ready-to-Use Set  
of Ophthalmic Gas

\*New Packaging Coming Soon

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE IFU SUPPLIED WITH THE PRODUCT

<https://www.bvimedical.com/customers-support/>



BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2024

1576806-02 IN